

3B Medical Lumin CPAP cleaner reviews

Lumin CPAP cleaner is the fastest and easiest way to clean and disinfect your CPAP mask, reservoir and other accessories without water, chemicals or harmful ozone. Lumin CPAP cleaner takes only five minutes to complete its clean cycle. It merely requires you to quickly wipe down your mask or water reservoir, place it into the Lumin's cleaning chamber, and turn it on with a press of a button. Just a short time later, the [Lumin CPAP cleaner](#) will chime, indicating that your equipment is clean and ready for immediate use.

How Does Lumin CPAP Mask Cleaner works?

3B Lumin CPAP cleaner and sanitizer device runs using a specialized high powered germicidal Ultra violet bulb. The drawer is made of polished aluminum which reflects the Ultraviolet light to improve protection. Ultra violet disinfection works on eliminating 99.99% bacteria, molds and germs and is well matched for keeping your CPAP mask and humidifier water chamber clean from germs that can make you sick.

Daily use of the Lumin CPAP cleaning machine will assist keep you healthy and compliant with your CPAP equipment. Because the 3B Lumin CPAP sanitizer device operates utilizing a type of high energy light, it is really efficient in disinfecting your CPAP accessories.

Introducing: Lumin™

The Best CPAP Cleaning System Available

The easiest automated CPAP sanitizer
Cleans and sanitizes CPAP mask and accessories

- No toxic liquids
- No heat or chemicals to clean your CPAP
- Automatic cycle time and LED indicator lights
- Simple to use and easy to clean
- The most advanced CPAP sanitizer

CPAP Mask Cleaner

CPAP Mask Cleaner

Lumin™

Lumin™

Lumin™

The Healthy Choice!
Ozone Free!
Kills 99% Harmful Germs

Learn More ▶

Main Features of Lumin CPAP cleaner

Safe and Easy to Utilize - The 5 minute ozone FREE cleaning cycle without any need for water or damaging chemicals.

Affordable - 3B Lumin CPAP and mask cleaning device budget friendly at less than other brands. A healthy CPAP experience ought to be budget friendly for everyone. It just cost \$249.00 & free shipping.

No Hazardous Ozone - UV light will disinfect approximately 99% of harmful mold and bacteria that can trigger infection and health problem. The usage of UV-C is eco-friendly, leaves no residue or poisonous gases or chemicals. UV-C systems are presently in usage to decontaminate ambulances, emergency situation service and other high touch places.

No Maintenance is needed - Lumin CPAP sanitizer needs no special care of maintenance and there are no parts that need to be replaced . The high powered UVC light bulb was developed specifically for the Lumin CPAP sanitizer and generates a high-powered disinfection wavelength of light and lasts longer than the life of your device, so no replacements or services are needed.

Multipurpose sanitizing machine - Lumin CPAP Cleaner is likewise a general-purpose disinfecting device and can disinfect any product put in the drawer. You can keep it on your bathroom counter and disinfect any essential items such as your toothbrush, child pacifiers and any toys.

Disinfect your CPAP equipment
in just 5 minutes!

AS EASY AS 1-2-3

- 1 Open the Lumin sliding tray
- 2 Place your CPAP mask and humidifier chamber on the Lumin tray
- 3 Close & push the button to start

UV kills 99% of bacteria. No harmful ozone! Only \$249.

Lumin™

BUY NOW

Lumin CPAP Cleaner - Frequently Asked Questions

See below the most frequently asked questions and answers below which help you to understand the product . For additional answers, feel free to call or email us.

Can I pay for Lumin using an HSA or FSA?

Yes, Lumin is eligible for purchase under both HSA and FSA.

How do I use the Lumin CPAP cleaner ?

Lumin CPAP sanitizer is a single button CPAP cleaner device. It really could not be any easier. First wipe the inside of your mask and cushion with a wipe to remove any physical residue. If cleaning your humidifier water chamber, open the chamber and place it in the drawer face up. If cleaning your mask, lay your mask in the drawer face up. Close the drawer and press the button. In 5 minutes the cycle is complete and your mask and water chamber are ready for use.

Is Lumin covered by Medicare or Insurance?

No, at this point in time, neither Medicare nor private insurers are covering cleaning devices for CPAP equipment.

How often do I need to clean my CPAP accessories?

You should maintain a proper CPAP cleaning schedule for your machine's accessories. 3B Medical recommends that you clean your CPAP accessories at least twice a week. A quick 5 minute cleaning cycle will ensure that your accessories are ready to use.

Does Lumin come with a warranty?

Yes, Lumin is backed by 3B Medical's comprehensive two year warranty.

Will Lumin damage my CPAP mask or accessories?

No. While long exposure times to UVC can weaken plastics and silicone, the Lumin cycle time is only 5 minutes. One year of daily use amounts to 30 hours of actual exposure time. In accelerated aging tests, 3B has validated no damage to plastics or silicone in over 150 hours of testing, roughly equivalent to 5 years of daily use.

Summary

Lumin CPAP cleaner is a wise investment for anybody who is looking for an easy to use and affordable CPAP cleaner . 3B Lumin is simple to operate and needs no maintenance. We recommend using Lumin CPAP sanitizer twice weekly on your mask, water chamber and tube.

Resources :

<https://www.youtube.com/watch?v=q4R9A47BhV0>

<https://sites.google.com/view/cpapcleaner/lumin-cpap-cleaner>